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---

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---

#### gDNA extraction from Arabidopsis for cloning

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## Rough gDNA extraction from Arabidopsis, Marchantia or Physcomitrium for quick genotyping

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**Solutions (both from sigmaaldrich.com):**

**溶液 (均购自 sigmaaldrich.com) :**

- (Extract-N-Amp Direct PCR Extraction Solution)
  - (Extract-N-Amp 直接 PCR 提取液)
- Dilution solution D5688-12mL (Dilution Solution, Plant Direct PCR Sigma-Aldrich)
  - 稀释液 D5688-12mL (植物直接 PCR 稀释液, Sigma-Aldrich)

Both solutions were originally intended to be used for the Extract-N-Amp Plant PCR kit (sigmaaldrich.com). However, a modified version of the manufacturer's protocol can be used to prep gDNA from different plant species for genotyping (up to 1000 bp fragments) and can be tested to be used for cloning (quality might not be good enough for difficult PCR amplifications or downstream cloning).

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**Protocol:**

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- Aliquot 15-30  $\mu$ L of extraction solution into a PCR tube.
  - 将 15-30  $\mu$ L 提取液分装到 PCR 管中。
- Transfer a small piece of plant material into the tube (for Arabidopsis: small leaf from seedlings or a piece of a leaf from older plants; for Marchantia or Physcomitrium: small piece of plant material).

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For genotyping (and subsequent sequencing): if a PCR volume of 50  $\mu$ L is used, 4  $\mu$ L of the gDNA should be recommended.

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---

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#### Synthetic genes:

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Currently, GenScript offers the best price for synthetic genes with the FLASH Gene service (link provided in original).

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Do not set-up an account directly; ask for an official recommendation from a colleague via the GenScript webpage, this will allow you to get some further discount for your quote.

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(If the prices change in the future, alternative companies could be: TwistBioscience or BioCat)

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Primers are ordered as “Single-stranded DNA” in tubes, generally at a concentration of 25 nmole DNA.

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### Reaction Setup for Phusion Polymerase:

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All components should be mixed and centrifuged prior to use. Assemble all reaction components on ice and quickly transfer the reactions to a thermocycler preheated to the denaturation temperature (98°C). Add the DNA Polymerase last to prevent any primer degradation caused by the 3'→5' exonuclease activity.

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### To determine the amount of template DNA:

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DNA	Amount (量)
Genomic (基因组)	50 ng–250 ng
plasmid or viral (质粒或病毒)	1 pg–10 ng

Note: do not use more DNA/plasmid template as indicated: too much template can consume all primers.

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### Thermocycling conditions for a routine PCR:

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### Thermocycling conditions (Touch-down example):

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### GeneJET PCR Purification Kit

### GeneJET PCR 纯化试剂盒

All purification steps should be carried out at room temperature. All centrifugations should be carried out in a table-top microcentrifuge at  $>12000 \times g$ .

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Note: 0.5 M EDTA is store-bought (do not attempt making this yourself). The 50 x TAE buffer is stored at room temperature. pH should not be adjusted normally.

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### DNA gel electrophoresis:

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- Weight out the appropriate mass of agarose into an Erlenmeyer flask. Agarose gels are prepared using a w/v percentage solution.
  - 称取适量的琼脂糖放入锥形瓶中。琼脂糖凝胶使用 w/v 百分比溶液制备。
  - The concentration of agarose in a gel will depend on the sizes of the DNA fragments to be separated, with most gels ranging between 0.5 %-2 %.
    - 凝胶中琼脂糖的浓度取决于要分离的 DNA 片段的大小, 大多数凝胶的浓度范围在 0.5%-2% 之间。
  - ▪ For several general approaches in our lab, a 1 % gel is appropriate
    - ▪ 对于我们实验室的几种常规方法, 1% 凝胶是合适的
  - ▪ For small DNA fragments (< 750 bp), try using a higher % (e.g. 2 %)
    - ▪ 对于小 DNA 片段 (< 750 bp), 尝试使用较高百分比 (例如 2%)
  - ▪ For separating large DNA fragments (> 5,000 bp), try using a lower % (e.g. 0.8 %)
    - ▪ 对于分离大 DNA 片段 (> 5,000 bp), 尝试使用较低百分比 (例如 0.8%)
- Add running buffer to the agarose (1 X TAE) inside the Erlenmeyer flask.
  - 向锥形瓶中的琼脂糖加入电泳缓冲液 (1 X TAE)。
  - ▪ Use around 30 mL for small gels
    - ▪ 小胶约用 30 mL
  - ▪ Use around 60 mL for medium gels
    - ▪ 中胶约用 60 mL
  - ▪ Use around 80 mL for large gels (two layers of combs)
    - ▪ 大胶 (双层梳子) 约用 80 mL
- Melt the agarose/buffer mixture in the microwave, make sure that the mixture is sufficiently dissolved (it can boil in the microwave).
  - 在微波炉中融化琼脂糖/缓冲液混合物, 确保混合物充分溶解 (它可以在微波炉中沸腾)。

- Pour the mixture into a prepared gel holder including the gel comb and add a droplet (1-2  $\mu\text{L}$ ) of Ethidium Bromide, mix thoroughly using the pipet tip.
  - 将混合物倒入准备好的插有梳子的制胶托盘中，并加入一滴 (1-2  $\mu\text{L}$ ) 溴化乙锭 (EtBr)，用移液枪头充分混合。
  - ▪ **BE CAREFUL** not to touch anything with your hands after working in the EtBr contaminated area! Always wear gloves when working with EtBr.
    - ▪ **小心**，在 EtBr 污染区工作后不要用手触摸任何东西！操作 EtBr 时始终佩戴手套。
  - ▪ Remove air bubbles attached to the comb wells.
    - ▪ 去除附着在梳孔上的气泡。
  - ▪ Make sure to have a tray underneath in case of spilling gel.
    - ▪ 确保下面有一个托盘，以防凝胶溢出。
- After solidification (30 min), place the gel into a gel chamber supplied with enough 1 x TAE running buffer.
  - 凝固后 (30 分钟)，将凝胶放入加有足够 1 x TAE 电泳缓冲液的电泳槽中。
  - ▪ Note: during warm summer days, gel solidification can take longer.
    - ▪ 注意：在炎热的夏季，凝胶凝固可能需要更长的时间。
- Gently remove the comb by lifting it upwards.
  - 通过向上提拉轻轻取下梳子。
- Load your samples.
  - 上样。
  - ▪ Samples contain loading dye to a 1 X end concentration.
    - ▪ 样品中含有终浓度为 1 X 的上样染料。
- Load a DNA ladder (e.g. 10 Kb-ladder: 4.5  $\mu\text{L}$ , at one or both sides of your sample).
  - 上载 DNA Ladder (例如 10 Kb-ladder: 4.5  $\mu\text{L}$ ，在样品的一侧或两侧)。
- Run the gel at manual settings (steady 120 V) for about 30-40 min or until your DNA is separated sufficiently.
  - 在手动设置下运行凝胶 (恒定 120 V) 约 30-40 分钟，或直到你的 DNA 充分分离。
- Visualize the gel at the Biorad GelDoc Gel Imaging System.
  - 在 Biorad GelDoc 凝胶成像系统中观察凝胶。

## RNA gel electrophoresis:

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Note: RNA is easily degraded by RNases which are omnipresent in the gel tank or the buffers. Also use filter tips to prevent aerosols being blown in your samples. Make sure all products are clean by rinsing them with RNase-AWAY. Wear gloves at all times.

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A full list of restriction enzymes available in the PLT lab can be found in “Documents - BIOSYST - Van de Poel Lab/PPL-Lists/Restriction enzymes PPL.xlsx”

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A typical restriction digest using ThermoFisherScientific restriction enzymes:

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Component	Amount
dH <sub>2</sub> O	Up to 20 $\mu$ L
10X Buffer	2 $\mu$ L (1X)
DNA	Up to 1 $\mu$ g
Restriction Enzyme (10 units)	1 $\mu$ L each

Mix gently and gently spin down for a few seconds. Note: the restriction enzyme should be the last added to the reaction and always kept on ice.

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- If you use several enzymes at once, the amount of enzymes should not exceed 1/10 of the total reaction volume.
  - 如果同时使用几种酶，酶的量不应超过总反应体积的 1/10。

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- If you use several enzymes at once, check for their efficiency in the diverse buffers (FD/B/G/O/R/Tango; Reaction Conditions for Restriction Enzymes | Thermo Fisher Scientific - BE).
  - 如果同时使用几种酶，请检查它们在不同缓冲液 (FD/B/G/O/R/Tango) 中的效率（参考 Thermo Fisher Scientific - BE 的限制性内切酶反应条件）。
- Incubation time can be 1-16 hours, typically 1-4 hours, but can be increased up to 16 hours with reduced units.
  - 孵育时间可以是 1-16 小时，通常为 1-4 小时，但如果减少酶单位，可以增加至 16 小时。
  - Longer incubation times are needed when buffer compatibility is not optimal.
    - 当缓冲液兼容性不是最佳时，需要更长的孵育时间。
- **Incubation temperature** is enzyme-dependent, typically at 37 °C.
  - 孵育温度取决于酶，通常为 37 °C。
- **Stopping the digestion** reaction might be done by heat inactivation (65 °C or 80 °C for 20 minutes, depending on the enzyme).
  - 终止消化反应可以通过热灭活（65 °C 或 80 °C 处理 20 分钟，具体取决于酶）来完成。

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Component	20 µL Reaction
T4 DNA Ligase Buffer (10X)*	2 µL
Vector DNA (50 ng)	Use NEBioCalculator to calculate molar ratios
Insert DNA	Use NEBioCalculator to calculate molar ratios
Nuclease-free water	to 20 µL
T4 DNA Ligase	1 µL

Note: You can easily scale down the reaction to 10 uL.

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T4 DNA Ligase Buffer should be thawed and resuspended at room temperature.

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- Gently mix the reaction by pipetting up and down and briefly spin down.
  - 通过上下吹打轻柔混合反应液并短暂离心。
- Incubate in a PCR machine:
  - 在 PCR 仪中孵育:
  - For cohesive (sticky) ends, incubate at 16 °C overnight or room temperature for 10 minutes.
    - ◦ 对于粘性末端，在 16 °C 过夜孵育或在室温下孵育 10 分钟。

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## Restriction-based cloning

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This is a classical method by combining the unique cut-sites of two different restriction enzymes to clone in a DNA fragment in the vector of choice. The RE cut-sites for a PCR product are typically added to the primers (extra base pairs after the template specific primer part).

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- Choose your backbone and insert so that they have compatible RE cut sites for restriction enzymes that allow your insert to be placed into the backbone in the proper orientation.
  - 选择你的骨架和插入片段，使它们具有兼容的限制性内切酶切割位点，从而允许你的插入片段以正确的方向放入骨架中。
- Note: take into account the orientation of the insert in your vector.
  - 注意：考虑载体中插入片段的方向。

- Depending on the restriction sites available in the backbone, add the corresponding sites to the start and end of the insert through PCR.
  - 根据骨架中可用的限制性位点，通过 PCR 将相应的位点添加到插入片段的起始和末端。
- Note: make sure these restriction sites are unique to the sites you need, avoid internal cut sites in both the vector and insert. In the case the RE can cut in your template, either use other RE or perform site directed mutagenesis.
  - 注意：确保这些限制性位点是你所需的唯一位点，避免载体和插入片段内部出现切割位点。如果酶在模板内部有切点，请使用其他酶或进行定点突变。
- Digest both backbone and insert with the correct RE(s) in separate reactions (see restriction digest).
  - 在单独的反应中用正确的限制性内切酶消化骨架和插入片段（见酶切消化）。
  - The restricted insert can be isolated by using a PCR purification kit (see PCR clean-up).
    - 酶切后的插入片段可以使用 PCR 纯化试剂盒分离（见 PCR 纯化）。
  - The restricted vector backbone will have to be gel extracted, therefore it is recommended to use 1 µg of DNA template for its digestion.
    - 酶切后的载体骨架必须进行胶回收，因此建议使用 1 µg DNA 模板进行消化。

Note: If you are going to use only one restriction enzyme for both ends, or enzymes that have compatible overhangs or no overhangs after digestion, you will need to use a phosphatase to prevent re-circularization of the backbone plasmid. Therefore, treat the vector with a phosphatase after RE digestion and prior to the ligation step (or prior to the gel purification step). CIP (calf alkaline phosphatase) or SAP (shrimp alkaline phosphatase) are commonly used phosphatases.

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## Gateway cloning

## Gateway 克隆

Official protocol for multiple Gateway cloning:

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The protocol is described for cloning a promotor and gene of interest using two entry vectors:

pDONRP4P1 for the promotor and pDONR221/pDONR207 for the coding sequence (CDS).

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## Primer design for cloning of promoter in pDONRP4P1:

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Note: the reverse primer should stop right before the ATG start site of the CDS.

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Note: in some cases, there are regulatory elements in the first intron of the gene, and in this case, you include the first intron in the sequence of the promoter. Make sure the start site is non-functional, so that the actual start site is in the fusion protein.

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## Primer design for cloning of coding sequences in pDONR221 or pDONR207:

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### Before starting:

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- For the BP and LR reactions, half the enzyme volume of the official protocol is sufficient (4  $\mu$ L instead of 8  $\mu$ L).
  - 对于 BP 和 LR 反应, 官方方案的一半酶体积就足够了 (4  $\mu$ L 代替 8  $\mu$ L)。
- Aliquot the BP and LR clonase enzyme mix and avoid freezing and thawing; prepare aliquots of about 2  $\mu$ L, and store at -80  $^{\circ}$ C as this increases the lifespan of the enzymes.
  - 分装 BP 和 LR 克隆酶混合物, 避免反复冻融; 准备约 2  $\mu$ L 的分装, 并储存在 -80  $^{\circ}$ C, 这可以增加酶的寿命。
- Treat the BP and LR clonase enzyme mix gently; do not vortex!
  - 轻柔对待 BP 和 LR 克隆酶混合物; **不要涡旋!**

### BP reaction (Gateway™ BP Clonase™ II Enzyme mix):

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- Incubate mixture at 25 °C for at least 1 h (preferably overnight) in a PCR machine.
  - 在 PCR 仪中于 25 °C 孵育混合物至少 1 小时 (最好过夜)。
  - ▪ Reaction mixture can be also kept at room temperature (instead of 25 °C), unless the room temperature is unusually high or low (winter/summer).
    - ▪ 反应混合物也可以保持在室温 (代替 25 °C) , 除非室温异常高或低 (冬天/夏天) 。
- Optional: Add 0.5 µL of Proteinase K solution to terminate the reaction; incubate for 10 min at 37°C.
  - \*works also without proteinase K treatment.
  - 可选: 加入 0.5 µL 蛋白酶 K 溶液以终止反应; 在 37°C 孵育 10 分钟。\*没有蛋白酶 K 处理也可以。

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For further details see: (Link in original)  
 更多详情请见: (原文中有链接)

### LR reaction (Gateway™ LR Clonase™ II Enzyme mix):

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Note: For multiple Gateway cloning:  
注意：对于多重 Gateway 克隆：

- LR reaction that combines multiple entry vectors into one destination vector needs to be done overnight.
  - 将多个入门载体组合到一个目的载体中的 LR 反应需要过夜进行。
- In this case, commercially produced competent *E. coli* cells are preferential.
  - 在这种情况下，商业生产的感受态大肠杆菌细胞是首选。

## GreenGate cloning

### GreenGate 克隆

GreenGate cloning is based on *Bsa*I mediated golden gate cloning, but adapted with a suite of plasmids for use in plants. The strategy consists of first making 6 entry plasmids (A – F modules, corresponding to promoter, N-terminal tag, CDS, C-terminal tag, terminator and a plant resistance gene) using *Bsa*I cloning.

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These entry modules are designed to give complementary overhangs to the neighboring plasmid in the scheme, so that they can only be ligated in the correct order A – F in the destination vector. 这些入门模块被设计为给方案中的相邻质粒提供互补的突出端，因此它们只能以正确的顺序 A – F 连接到目的载体中。

Then, these 6 plasmids and the destination plasmid are taken to a single tube Greengate reaction, where they are digested with *Bsa*I and ligated with T4 DNA ligase. The standard GreenGate protocol results in a single expression cassette, but it is also possible to have multiple expression cassettes.

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The benefits of Greengate are that it's modular, scar-free and easy to make single transcriptional units. The downside is that it's not ideal if the regions you are coding have internal *Bsa*I sites, as ideally this would require mutagenesis to remove them before use.

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We have the original GreenGate kit (Addgene kit 1000000036; <https://www.addgene.org/kits/lohmann-greengate/>; Lampropoulos A, Sutikovic Z, Wenzl C, Maegele I, Lohmann JU, et al. (2013) GreenGate - A Novel, Versatile, and Efficient Cloning System for Plant Transgenesis. PLOS ONE 8(12): e83043.

<https://doi.org/10.1371/journal.pone.0083043>). The full list of plasmids is available from Addgene or the shared –80 excel sheet. We also have added GFP in the B module.

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**GreenGate vector design and layout.** The GreenGate system six modules represent the plant promoter (1), an N-terminal tag (2), the coding sequence of the gene of interest (3), a C-terminal tag (4), the plant terminator (5) and the plant resistance cassette (6).

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### • GreenGate entry vector (入门载体)

- The GreenGate cloning system uses six different types of pUC19 based entry vectors into which the individual elements are inserted.
  - GreenGate 克隆系统使用六种不同类型的基于 pUC19 的入门载体，单个元件被插入其中。
- The multiple cloning site of pUC19 has been replaced by two *Bsal* recognition sites (magenta scissors), the respective overhangs for each module type, and a counter-selectable *ccdB* gene. DNA fragments can be cloned via the specific overhangs, via the BamHI and KpnI sites or via A-overhangs after XcmI digestion.
  - pUC19 的多克隆位点已被两个 *Bsal* 识别位点（品红色剪刀）、每种模块类型各自的突出端以及一个反向筛选的 *ccdB* 基因所取代。DNA 片段可以通过特定的突出端、BamHI 和 KpnI 位点或 XcmI 消化后的 A 突出端进行克隆。

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### • GreenGate destination vector (目的载体)

- The GreenGate cloning system uses a pGreen-IIS based destination vector.
  - GreenGate 克隆系统使用基于 pGreen-IIS 的目的载体。
- The empty destination vector consists of a counter-selectable *ccdB*-cassette between the LB and RB sequences of pGreen-IIS, flanked by *Bsal* sites, with overhangs A and G.
  - 空的目的载体由 pGreen-IIS 的 LB 和 RB 序列之间的反向筛选 *ccdB* 盒组成，两侧是 *Bsal* 位点，具有 A 和 G 突出端。
- Note: The pSa origin of replication (ori A. tum.) requires the presence of the helper plasmid pSOUP in agrobacteria.
  - 注意：pSa 复制起始位点 (ori A. tum.) 需要农杆菌中辅助质粒 pSOUP 的存在。

### Making entry modules:

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PCR is used to amplify the region of interest and add appropriate *Bsal* sites for cloning. The following sequences are added in front of the gene specific region of the F primer, and in front of the reverse complement of the gene specific region for the R primer.

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**F:** 5' AACAGGTCTCANNNN (nn) - Gene specific region -3'

**R:** 5'-AACAGGTCTCANNNN - Gene specific region reverse complement 3'

- GGTCTC is the *Bsal* recognition site.
  - GGTCTC 是 *Bsal* 识别位点。

- AACAA was added because the enzyme does not cut if the restriction site is at the extreme ends of PCR products.
  - 添加 AACAA 是因为如果限制性位点位于 PCR 产物的极末端，酶不会切割（保护碱基）。
- NNNN represents the module specific overhang with the coding triplet being underlined.
  - NNNN 代表模块特定的突出端，编码三联体带有下划线。
- and 2 nucleotides (nn) are needed in case of the coding sequence and C-tag modules to bring the modules into frame.
  - 对于编码序列和 C 标签模块，需要 2 个核苷酸 (nn) 使模块进入读码框。

PCR is conducted using standard protocols for high fidelity polymerase, and PCR product purified by gel extraction.

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Both PCR product and the corresponding entry module is digested with *BsaI*, and purified (PCR purification kit for PCR product, preferably gel extraction for entry module).

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Then the digested entry PCR product and the corresponding entry vectors are ligated.

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The resulting entry vectors are then taken to the GreenGate reaction.

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Cycle at 37 C for 2 mins and 16 C for 2 mins for 30 cycles, then 80 C for 5 mins. Transform 5 ul of the reaction into appropriate cells.

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## Multiple transcriptional units

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#### Two transcriptional units

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For two transcriptional units, two intermediate vectors, pGGM000 and pGGN000, and two adapter molecules are used.

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The final construct is created in two steps: First, the two transcriptional units are assembled separately in two parallel reactions into the intermediate vectors.

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The first of these two supermodules consists of a plant promoter, an N-terminal tag, a coding sequence, a C-terminal tag and a plant terminator.

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Furthermore, a pre-cloned adapter module with F and H overhangs is added to the reaction.

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The intermediate vector has *Bsa*I recognition sites remaining in the vector backbone with matching A and H overhangs and confers resistance to kanamycin.

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The elements are combined in a GreenGate reaction but because two *Bsa*I sites remain in the final construct, an additional ligation reaction is performed before transformation.

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The second supermodule is built similarly, but the adapter molecule has H and A overhangs and at the 39-end a plant resistance cassette is added.

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Thus, the resulting plasmid will have H and G overhangs.

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The two supermodules are then combined with the destination vector in a second step in a standard GreenGate reaction yielding the final construct.

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## More than two transcriptional units

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An additional module that is a synthetic DNA duplex with unpaired ends, allows ligation to H and G overhangs during a GreenGate reaction to make it possible to have two transcriptional units.

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It contains two internal *Bsa*I recognition sites, which cannot be cut by *Bsa*I due to overlapping cytosine methylation.

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Thus, the element will behave as an ordinary H to G module in the first GreenGate reaction and replace the resistance cassette.

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After transformation and miniprepping, this methylation will be lost as there are no dcm sites present in the oligo duplex.

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In a second round of GG reaction, the two internal *Bsa*I recognition sites are now targets for *Bsa*I, giving rise to A and G overhangs, simulating an empty destination vector.

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Therefore, in each subsequent GreenGate round either the uncleavable oligo duplex can be used to prepare the plasmid for the insertion of another expression cassette or the plasmid can be finalized by adding a plant resistance cassette to the reaction.

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In this approach, the transcriptional units are assembled serially, so it requires one additional step for each expression cassette to be added.

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The FHadapter used along with the oligo duplex separates the individual expression units in the final construct and thus can be freely designed to act as a spacer of arbitrary length and sequence minimizing the risk of mutual interference of the expression behavior.

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Protocol Samuel Miller Lab

Samuel Miller 实验室方案

Gibson assembly primers are broken down in two parts: primer sequence and overlap sequence. The primer sequence should be designed using traditional characteristics in mind (i.e.  $T_m$  values, G/C ratio, and G/C anchors). The overlap sequence needs to have between 20 – 150 bp homology to insert or vector. A good recommended start is a primer of 60 bp in total length: 30 bp of vector and 30 bp of insert.

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## Gibson 反应

### Reaction Stock Preparation

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Prepare 5x ISO buffer: 6 mL of this buffer can be prepared by combining the following:

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Aliquot in 100  $\mu$ L aliquots and store at -20 °C.

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### Prepare master mixture. Combine the following:

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Aliquot in 15  $\mu$ L aliquots and store at -20 °C.

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This assembly mixture can be stored at -20 °C for at least one year.

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### Reaction Procedure

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- Thaw a 15  $\mu$ L assembly mixture aliquot and keep on ice until ready to be used.

- 解冻一份 15  $\mu$ L 组装混合物，并在使用前保持在冰上。
- Add 5  $\mu$ L of PCR-amplified vector and insert DNA to be assembled to the master mixture.
  - 向母液中加入 5  $\mu$ L PCR 扩增的载体和要组装的插入 DNA。
  - Note: The DNA should be in equimolar amounts.
    - 注意：DNA 应为等摩尔量。
  - Total DNA to be assembled should range between 20- 200 ng.
    - 组装的总 DNA 量应在 20-200 ng 之间。
  - For large DNA segments (> 150 kb): increased DNA amount should be used (e.g. 250 ng of each 150 kb DNA segment).
    - 对于大 DNA 片段 (> 150 kb): 应使用增加的 DNA 量 (例如每个 150 kb DNA 片段使用 250 ng) 。
- Thoroughly mix gently.
  - 轻柔地充分混合。
- Incubate at 50 °C for 15 to 60 min.
  - 在 50 °C 孵育 15 到 60 分钟。
- If cloning is desired, transform 1 – 5  $\mu$ L to competent E. coli cell.
  - 如果需要克隆，将 1 – 5  $\mu$ L 转化到大肠杆菌感受态细胞中。
  - Note: Dialyze prior to electroporation. Heat-shock transformation can be done right so.
    - 注意：电穿孔前需要透析 (去除盐分) 。热激转化可以直接进行。

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A major limitation (especially in high GC-rich templates such as in *P. aeruginosa*) to Gibson assembly is that the termini of the DNA sequence fragments to be assembled, should not have stable single stranded DNA secondary structure, such as a hairpin or a stem loop, or repeated sequences, as this would directly compete with the required single-stranded annealing/priming of neighboring assembly fragments.

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So make sure to run your sequences in Scitools to look for such nuisance areas (Delta G for hairpins should be more than -10 and with a lower annealing temperature than that of primer:primer duplex). 因此，请确保在 Scitools 中运行您的序列以查找此类麻烦区域 (发夹的 Delta G 应大于 -10，且退火温度低于引物:引物双链体的退火温度) 。

## Site directed mutagenesis (SDM)

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#### Blunt end ligation-based SDM methods:

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SDM using blunt end ligation results in a blunt end PCR product, which is circularized by ligation. First, primers are designed to incorporate the desired mutation (substitution, deletion or insertion). The plasmid to mutagenize is used as a template for PCR, and the resulting blunt end PCR product is ligated.

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Because the PCR product is missing the 5' phosphate, this needs to be added using a kinase. To avoid returning the original vector used as a template in the PCR, the sample should be digested using *DpnI*, which only cleaves methylated DNA and not the PCR product. Finally, the blunt end PCR product should be circularized using a ligase.

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Component	Amount
DNA	Up to 1 µg
10X FastDigest buffer	2 µL
<i>DpnI</i>	1 µL
Nuclease free H <sub>2</sub> O	To 20 µL

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Component	Amount
DNA	Up to 300 pmol of 5' termini*
T4 PNK Reaction Buffer 10X	5 µL
ATP 10 mM	5 µL
T4 PNK	1 µL 10 units
Nuclease free H <sub>2</sub> O	To 50 µL

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The T4 PNK buffer can be substituted for T4 DNA ligase buffer. The 1X T4 DNA ligase buffer contains 1 mM ATP, so ATP does not need to be added. In this case, the product does not need to be cleaned before ligation, but supplemental ATP should be added to the kinase reaction to a final concentration of 1 mM.

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Component	Amount
DNA	1 µL of 50 ng
T4 Ligase Buffer 10X	2 µL
T4 Ligase	1 µL
Nuclease free H2O	To 20 µL

- o Incubate at 16 °C overnight or at room temperature for 10 minutes.
  - o o 在 16 °C 过夜孵育或在室温下孵育 10 分钟。
  - o ▪ The combination of 10 min at room temperature and 16 °C overnight is also possible.
    - ▪ 室温 10 分钟和 16 °C 过夜的组合也是可能的。
- o Heat inactivate at 65°C for 10 mins.
  - o o 在 65°C 下热灭活 10 分钟。
- o Transform 1-5 µL into appropriate E. coli strain.
  - o o 将 1-5 µL 转化到适当的大肠杆菌菌株中。

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Component	Amount
PCR product	1 $\mu$ L
2X KLD reaction buffer	5 $\mu$ L
10X KLD enzyme mix	1 $\mu$ L
Nuclease free H <sub>2</sub> O	3 $\mu$ L

- ◦ Mix well by pipetting up and down and incubate at room temperature for 5 minutes.
  - ◦ 通过上下吹打充分混合并在室温下孵育 5 分钟。
- ◦ Transform 5  $\mu$ L of KLD reaction to 50  $\mu$ L appropriate competent cells.
  - ◦ 将 5  $\mu$ L KLD 反应液转化到 50  $\mu$ L 适当的感受态细胞中。

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Gibson based SDM also requires designing primers incorporating the desired mutation, but the resulting PCR product should have complementary overlapping ends of 15-30 bp. Then an exonuclease chews back the 5' end, complementary regions anneal, DNA polymerase fills in the gaps within each annealed fragment, and ligase seals the nicks. The benefit of this approach is that it can also be used to assemble multiple fragments in one go.

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### Recommended Amount of Fragments Used for Assembly

Component	2-3 Fragment Assembly*	4-6 Fragment Assembly**
Recommended DNA Molar Ratio	vector:insert = 1:2	vector:insert = 1:1
Total Amount of Fragments	0.03-0.2 pmols* (X $\mu$ L)	0.2-0.5 pmols** (X $\mu$ L)
NEBuilder HiFi DNA Assembly Master Mix	10 $\mu$ L	10 $\mu$ L
Nuclease free H <sub>2</sub> O	10-X $\mu$ L	10-X $\mu$ L
<b>Total Volume</b>	<b>20 <math>\mu</math>L ††</b>	<b>20 <math>\mu</math>L ††</b>

- - o Optimized cloning efficiency is 50-100 ng of vector with 2-fold excess of each insert. Use 5-fold molar excess of any insert(s) less than 200 bp. Total volume of unpurified PCR fragments in the assembly reaction should not exceed 20%. To achieve optimal assembly efficiency, design 15-20 bp overlap regions between each fragment.
    - 优化的克隆效率是用 50-100 ng 载体和 2 倍过量的每个插入片段。对于小于 200 bp 的任何插入片段，使用 5 倍摩尔过量。组装反应中未纯化 PCR 片段的总量不应超过 20%。为获得最佳组装效率，在每个片段之间设计 15-20 bp 的重叠区域。
  - o \*\* To achieve optimal assembly efficiency, design 20-30 bp overlap regions between each fragment with equimolarity of all fragments (suggested: 0.05 pmol each).
    - o \*\* 为获得最佳组装效率，在每个片段之间设计 20-30 bp 的重叠区域，所有片段等摩尔（建议：每个 0.05 pmol）。
  - o † Control reagents are provided for 5 experiments.

- † 提供的对照试剂可供 5 次实验使用。
- †† If greater numbers of fragments are assembled, increase the volume of the reaction, and use additional NEBuilder HiFi DNA Assembly Master Mix.
  - †† 如果组装更多数量的片段，请增加反应体积，并使用额外的 NEBuilder HiFi DNA 组装母液。
- ○ Incubate samples in a thermocycler at 50 °C for 15 minutes (when 2 or 3 fragments are being assembled) or 60 minutes (when 4–6 fragments are being assembled).
  - ○ 在热循环仪中 50 °C 孵育样品 15 分钟（组装 2 或 3 个片段时）或 60 分钟（组装 4-6 个片段时）。
- ○ Store samples on ice or at –20 °C for subsequent transformation.
  - ○ 将样品储存在冰上或 –20 °C 以备后续转化。
- ○ Transform 2-5 uL of assembled product to appropriate cell strain.
  - ○ 将 2-5 uL 组装产物转化到适当的细胞菌株中。

### Restriction enzyme-based SDM methods:

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Restriction enzyme based SDM makes use of type II restriction enzymes (we have *Bsal*, *AarI* and *BbsI*). These enzymes cut outside of their recognition site, meaning cloning can be scar-free.

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Primers can be designed to introduce the appropriate mutation with internal restriction sites. This can be designed so that the recognition sites are cut out, leaving complementary sticky ends which are ligated. The benefit of this method is that multiple products can be ligated together, but it is best if there are no internal sites for the restriction enzymes which are used.

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### Type II restriction enzyme mutagenesis

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Component	Amount
dH <sub>2</sub> O	Up to 20 μL
10X Buffer*	2 μL (1X)
DNA	Up to 1 μg
Restriction Enzyme (10 units)	1 μL each

- \*If using Thermo *AarI*, be sure to add the 50X oligonucleotide mix as well as the buffer.
  - \*如果使用 Thermo *AarI*, 请务必添加 50X 寡核苷酸混合物以及缓冲液。
- ◦ PCR clean PCR product
  - ◦ 纯化 PCR 产物
- ◦ Ligate using T4 DNA ligase
  - ◦ 使用 T4 DNA 连接酶进行连接

Component	Amount
DNA	1 μL of 50 ng

Component	Amount
T4 Ligase Buffer 10X	2 $\mu$ L
T4 Ligase	1 $\mu$ L
Nuclease free H <sub>2</sub> O	To 20 $\mu$ L

- o Transform 1  $\mu$ L of ligation to appropriate cells.
  - o 将 1  $\mu$ L 连接产物转化到适当的细胞中。

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### Competent cells available in the lab:

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## Agrobacterium tumefaciens transformation

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#### Transformation of electrocompetent *A. tumefaciens* cells using electroporation:

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- Carefully pipet *A. tumefaciens* and plasmid DNA mixture up and down for 1-2 times.
  - 小心上下吹打 1-2 次混合农杆菌和质粒 DNA。
- Transfer cell-plasmid-mixture into pre-cooled (on ice) electroporation cuvettes.
  - 将细胞-质粒混合物转移到预冷（冰上）的电击杯中。
- Wipe ice/any residual water from the cuvettes before transferring it into the electroporator.

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- PH to 7.0 with NaOH (用 NaOH 调节 pH 至 7.0)
- Autoclave (高压灭菌)

#### **50 mL 1x YT broth:**

#### **50 mL 1x YT 肉汤:**

- 0.4 g Bacto tryptone
- 0.25 g Bacto yeast extract
- 0.125 g NaCl
- Water up to 50 mL
- pH to 7.0 with NaOH
- Autoclave

#### **TS buffer (TSB; prepare fresh on the day, when the competent cells are being made):**

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#### **For all steps, work in a laminar flow to keep the cells sterile**

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## Day 3:

### From here on, keep the cultures/cells always on ice

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For all the steps, work in a laminar flow to keep everything sterile.

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## Colony PCR

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Colony PCRs can be used as an initial test to confirm success of the cloning approach before inoculating the transformants in liquid culture. For this, select a region for not more than 1 kb that you want to amplify. The primers that you use for amplifying the region of interest should be insert- and vector-specific (e.g. forward primer binds within the insert, reverse primer binds in the vector backbone).

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Depending on the number of colonies you want to test, prepare a master mix with the ThermoScientific DreamTaq Green PCR Master Mix (2X).

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It is not necessary to prepare reaction volumes of up to 50  $\mu\text{L}$  (according to the manufacturer's instructions), one quarter of the reaction volume (12.5  $\mu\text{L}$ ) is sufficient.

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**Always prepare a master mix for 1 reaction in excess.**

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**Adding material of the bacterial colonies into the PCR vials containing the PCR mixture:**

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- Aliquot the PCR mixture into PCR vials and transfer those into the laminar flow.
    - 将 PCR 混合物分装到 PCR 管中，并将其转移到层流罩中。
  - Touch the bacterial colony with a sterile tip and transfer it into the PCR vials; gently pipette up and down to resuspend the bacteria in the PCR mixture.
    - 用无菌枪头接触细菌菌落并将其转移到 PCR 管中；轻轻上下吹打以将细菌重悬在 PCR 混合物中。
    - It is not necessary to visually see bacteria on the pipette tip, it is sufficient to gently touch the colony.
      - 不需要在移液枪头上肉眼看到细菌，只需轻轻接触菌落即可。
- 💡 注释:** 菌量过多会抑制 PCR 反应。
- **Tip: Streak the colonies you want to test on a “Master” plate before adding it to the PCR mixture; then, you have more material for later inoculations and can organize them more carefully.**
    - **提示: 在将其加入 PCR 混合物之前, 将你要测试的菌落在“主”平板上划线 (备份); 这样, 你就有了更多材料供以后的接种使用, 并且可以更仔细地组织它们。**
  - Spin down the PCR reaction mixture and start the PCR in a thermo cycler.
    - 离心 PCR 反应混合物并在热循环仪中开始 PCR。

**PCR settings:**

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- Perform the PCR according to the manufacturer's instructions, to the length of the fragment you want to amplify and the annealing temperature of the primers.

- 根据制造商的说明、你要扩增的片段长度和引物的退火温度进行 PCR。
- 30 cycles are enough for colony PCRs; otherwise, you may proliferate artificial products.
  - 30 个循环对于菌落 PCR 来说足够了；否则，你可能会扩增出人造产物（非特异性产物）。
- **E. coli colony PCRs: for initial denaturation step at 98 °C for 3 min.**
  - **大肠杆菌菌落 PCR: 初始变性步骤 98 °C 3 分钟。**
- **A. tumefaciens colony PCRs: for initial denaturation step at 98 °C for 10 min.**
  - **农杆菌菌落 PCR: 初始变性步骤 98 °C 10 分钟。**

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After performing the PCR reaction in a thermo cycler, run the PCR in on a 1 % agarose gel together with a 1 kB DNA ladder. Afterwards, visualize the DNA bands with UV light. If the DNA bands are in accordance with the expected fragment size, you can inoculate the bacterial colony into LB medium with the respective antibiotic for further downstream processing (e.g. plasmid prep).

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## Sequencing (Sanger)

### 测序 (Sanger 法)

The results of our cloning approaches can be verified by Sanger Sequencing. We are sending our samples to LGC Genomics (Link in original). Results can be analyzed using Geneious.

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#### Step 1: Prepare your samples in 1.5 mL reaction tubes

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## Step 2: Place your order at LGC

### 第 2 步: 在 LGC 下单 (Link in original)

- At "sequencing ordering", you will find below "single sample" the option for "Ready2 / Flexi Run".
  - 在 "测序订购" 处, 你会在 "单一样品" 下方找到 "Ready2 / Flexi Run" 选项。
- Enter your sample label ID, sample name and the service you want to use (Ready2 / Flexi Run).
  - 输入你的样品标签 ID、样品名称和你想要使用的服务 (Ready2 / Flexi Run)。
  - In case you chose Flexi Run, indicate which primer should be used by LGC.
    - 如果你选择 Flexi Run, 请注明 LGC 应使用哪种引物。
- Continue with filling in the information for the pick-up box and your e-mail address on the next site, and confirm the order.
  - 继续在下一个站点填写取件箱信息和你的电子邮件地址, 并确认订单。
- Schedule the pick-up with DHL until 12:00 via (Link provided).
  - 通过 (提供的链接) 在 12:00 之前安排 DHL 取件。
- Results will be available in about 2-3 working days (notification via e-mail) and can be downloaded at "Data Management – Order&Download Results".
  - 结果将在大约 2-3 个工作日内提供 (通过电子邮件通知), 并可在 "数据管理 – 订单和下载结果" 处下载。

## Step 3: Analyze your sequences with Geneious

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---

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---

- Prepare stock solutions of each antibiotic according to the desired concentration.
  - 根据所需浓度制备每种抗生素的储备液。
- Filter-sterilize in the laminar flow and aliquot into 1.5 mL tubes.

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**For solid plates: Add antibiotics to the medium only, if it is hand-warm or cooler!**

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